


## Optimum Wound Approximation & Healing

**MIRUS™** Disposable Skin Stapler with 316L stainless steel staples has applications in abdominal, gynecological, orthopedic, and thoracic surgeries for skin closure. Also, it could be used for variety of skin closures. After wound approximation and healing, these stainless steel staples are then removed by using a Skin Staple Extractor.

Features:	Benefits:
Ergonomic design	Easy handling for desired result
Optimum leg length	Clasp thicker tissues
Transparent housing with Indicator	Enhance staple count visibility
Adequate staple crown length	Secures wider wounds



**Staple material:** 316L Stainless steel  
**Sterilisation:** Ethylene oxide  
**Packaging:** 6 Pcs/Box

**Code:** MSS35  
**Size:** 35  
**Staples:** Wide

Larger diameter wire, greater width and extra staple height lead to secure wound closure